

## MYCOTOXINS

Rapid Liquid Chromatographic Determination of Aflatoxins  $M_1$  and  $M_2$  in Artificially Contaminated Fluid Milks: Collaborative StudyROBERT D. STUBBLEFIELD and WILLIAM F. KWOLEK<sup>1</sup>*U.S. Department of Agriculture, Agricultural Research Service, Northern Regional Research Center, Peoria, IL 61604*

Collaborators: R. Beaver; M. Billotte; J. Ferguson-Foos; J. I. Greer; K. D. Johnson; R. Lyon; E. H. Marth; J. I. Pitt; M. E. Rhodes; J. L. Richard; B. Roberts; G. Rusul; Y. Saito; E. A. Sizoo; G. Sutton; N. F. Tobin; M. Trucksess; H. P. van Egmond; J. D. Warren; N. Webb; D. Wilson; A. E. Yousef

An international collaborative study involving 14 collaborators from 5 different countries was conducted to test a rapid liquid chromatographic (LC) method for detecting aflatoxins  $M_1$  and  $M_2$  in fluid milk. Each collaborator prepared artificially contaminated milk samples (0.078–1.31 ng  $M_1$ /mL and 0.030–0.13 ng  $M_2$ /mL) by adding solutions containing various concentrations of aflatoxins  $M_1$  and  $M_2$  to fresh milk. Recoveries ranged from 85.2 to 102.5% (av. 93.7%) for aflatoxin  $M_1$  and from 99.5 to 126.7% (av. 109.8%) for aflatoxin  $M_2$ . Coefficients of variation averaged 21.4% ( $M_1$ ) and 35.9% ( $M_2$ ). An analysis of variance was calculated from combined data to determine variance components. The within-laboratory variations ( $S_o$ ) (repeatability) were 27.9% ( $M_1$ ) and 23.9% ( $M_2$ ), and the among-laboratory variations ( $S_b$ ) (reproducibility) were 44.5% ( $M_1$ ) and 64.7% ( $M_2$ ). No visual differences were determined between normal or reverse phase LC for contaminated samples; however, there were an insufficient number of collaborators using normal phase to give meaningful separate statistical data. For 26 observations of uncontaminated milk, 3 false  $M_1$  positives were reported for normal phase LC determinations and 2 false  $M_1$  positives were reported for reverse phase LC determinations. Three normal phase and 11 reverse phase false  $M_2$  positives were reported for 104 observations in uncontaminated milk. The reverse phase LC method for determination of aflatoxins  $M_1$  and  $M_2$  in fluid milk has been adopted official first action.

In 1979, an international collaborative study (1) was conducted to test the Stubblefield method (2) for determining aflatoxin  $M_1$  in dairy products and the van Egmond et al. method (3) for thin layer chromatographic (TLC) confirmation of  $M_1$  identity. Subsequently, both methods were adopted by AOAC and the Commission of Food Contaminants of the International Union of Pure and Applied Chemistry (IUPAC) (4). The analytical method (2) uses TLC for the quantitation step. Several methods that incorporate liquid chromatography (LC) for determining  $M_1$  in milk have been published (5–11). Some of these methods utilize disposable silica gel and/or bonded silica gel ( $C_{18}$ ) cartridges to reduce analysis time. Minimum detection limits for the LC methods range from 0.01 to 0.10 ppb aflatoxin  $M_1$  and can be adjusted readily within this range as required by the guidelines or tolerances of state, federal, or international agencies.

The Associate Referee evaluated each method to determine which one(s) provided: (1) interference-free extracts; and (2) satisfactory recoveries (>80%) of added  $M_1$ . From the data

and comments of other scientists (personal communications), the Foos and Warren method (11) was selected for international collaborative study. This method, as published, was developed for normal phase LC. It is desirable that the official AOAC method be applicable to reverse phase LC, also. The published procedures for reverse phase LC of aflatoxin  $M_1$  (6–9) were investigated, and a modification of the Beebe and Takahashi method (6) was selected. This method measures the fluorescent derivative formed by treating  $M_1$  with trifluoroacetic acid (TFA). In aqueous mobile phase solvents, the derivative is more fluorescent than  $M_1$  itself.

In 1983, an international collaborative study was initiated; however, incomplete reaction of standard aflatoxin  $M_1$  and TFA occurred, and the study was halted. These problems were essentially eliminated by forming the standard derivative in a silylated glass vial to prevent irreversible adsorption of the  $M_1$ – $M_2$  standards to the glass walls (future publication). Adsorption causes incomplete reaction between  $M_1$  and TFA. Further research by the Associate Referee has shown that an increase in the reaction temperature (40°C) is necessary to achieve complete derivatization.

All attempts to find a suitable preservative for naturally contaminated fluid milk were unsuccessful. Chemicals were found which preserved the milk; however, the  $M_1$  in the milk was degraded. Consequently, collaborators were asked to furnish uncontaminated milk and to prepare artificially contaminated milk samples with sealed acetonitrile solutions of aflatoxins  $M_1$  and  $M_2$ . The report of the data submitted by 14 collaborators from 5 different countries is presented here.

## Collaborative Study

Aflatoxin  $M_1$ – $M_2$  Standard Solutions

Crystalline aflatoxins  $M_1$  and  $M_2$  were used to prepare stock solutions of each aflatoxin (266.07  $\mu$ g  $M_1$ /mL and 172.02  $\mu$ g  $M_2$ /mL, in acetonitrile). Aflatoxin concentrations in the stock solutions were determined according to 26.004–26.011 (12), using extinction coefficients of 19 850 and 21 400 for  $M_1$  and  $M_2$ , respectively, in acetonitrile. Purity criteria for crystalline  $M_1$  and  $M_2$  are given by Stubblefield et al. (13, 14). The stock solutions were used to prepare a standard solution for LC, containing 0.50  $\mu$ g  $M_1$  and 0.10  $\mu$ g  $M_2$ /mL in acetonitrile–benzene (1 + 9).

## Preparation of Samples

All samples in the study were artificially contaminated fluid milks. Seven aflatoxin  $M_1$ – $M_2$  spiking solutions (in acetonitrile) were prepared from aliquots of  $M_1$  and  $M_2$  stock or diluted stock solutions in separate 100 mL volumetric flasks as follows: Samples 1/8, and practice, 29.1  $\mu$ L stock  $M_1$  with 37.1  $\mu$ L 10-fold dilution  $M_2$ ; samples 2/9, 58.2  $\mu$ L 10-fold dilution stock  $M_1$ ; samples 3/10, 29.2  $\mu$ L 10-fold dilution stock

Submitted for publication October 1, 1985.

This report of the Associate Referee, R. D. Stubblefield, was presented at the AOAC 10th Annual Spring Training Workshop, April 9–11, 1985, at Dallas, TX.

The recommendation of the Associate Referee was approved by the General Referee and the Committee on Foods I and was adopted by the Association. See the General Referee and Committee reports, *J. Assoc. Off. Anal. Chem.* (1986) 69, March issue.

The mention of firm names or trade products does not imply that they are endorsed or recommended by the U.S. Department of Agriculture over other firms or similar products not mentioned.

<sup>1</sup>Retired.

Table 1. Collaborative results for determination of aflatoxins M<sub>1</sub> and M<sub>2</sub> (ng/mL) in artificially contaminated fluid milk (blind duplicate pairs)\*

Coll.	Method <sup>b</sup>	Sample 1		Sample 8		Sample 2		Sample 9		Sample 3		Sample 10	
		M <sub>1</sub>	M <sub>2</sub>	M <sub>1</sub>	M <sub>2</sub>	M <sub>1</sub>	M <sub>2</sub>	M <sub>1</sub>	M <sub>2</sub>	M <sub>1</sub>	M <sub>2</sub>	M <sub>1</sub>	M <sub>2</sub>
1 <sup>c</sup>	R	0.32	0.01	0.20	0.03	0.05	0	0.05	0	0.05	0	0.03	0
4	R	0.129 <sup>d</sup>	0.010	0.287	0.023	0.112	0	0.132	0	0.073	0	0.034	0
6	R	0.85	0.07	1.07	0.10	0.20	0.02	0.16	0	0.06	0	0.06	0
8	N	0.80	0.10	0.81	0.09	0.12	0.017	0.10	0	0.06	0	0.06	0
18	R	0.72	0.05	0.68	0.06	0.22	0.01	0.10	0	0.07	0	0.08	0
19	R	0.76	0.16	0.18	0	0.074	0.18	0.19	0	0.048	0.17	0.081	0.082
20	N	0.85	0.08	0.83	0.08	0.17	0	0.21	tr <sup>e</sup>	0.07	0	0.11	0
21	N	0.7	0	0.7	0	0.3	0	0.3	0	0.2 <sup>d</sup>	0	0.2 <sup>d</sup>	0
22	R	0.81	0.1	0.73	0.08	0.20	0	0.20	0	0.11	0	0.12	0
23	R	0.835	0.084	0.925	0.086	0.165	0	0.085	0	0.095	0.02	0.076	0.05
25	R	0.57	0.077	0.71	0.055	0.15	0	0.13	0	0.063	0	0.046	0
26	R	0.77	0.06	0.73	0.08	0.18	0	0.18	0	0.09	0	0.10	0
27	R	0.49	0.03	0.30	0.05	0.07	0	0.10	0	0.02	0	0.02	0
33	R	0.75	0.09	0.70	0.08	0.25	0	0.15	0	0.20 <sup>d</sup>	0	0.20 <sup>d</sup>	0

Coll.	Method <sup>b</sup>	Sample 5		Sample 12		Sample 6		Sample 13		Sample 7		Sample 14	
		M <sub>1</sub>	M <sub>2</sub>	M <sub>1</sub>	M <sub>2</sub>	M <sub>1</sub>	M <sub>2</sub>	M <sub>1</sub>	M <sub>2</sub>	M <sub>1</sub>	M <sub>2</sub>	M <sub>1</sub>	M <sub>2</sub>
1 <sup>c</sup>	R	0.08	<0.01	0	0	0.30	0.05	0.26	0.03	0.06	0	0.02	0
4	R	0.360	0.026	0.183	0.015	0.815	0.075	0.512	0.055	0.069	0	0.044	0
6	R	0.54	0.05	0.70	0.04	1.59	0.145	1.71	0.18	0.17	0.01	0.14	0
8	N	0.44	tr <sup>e</sup>	0.51	0.05	1.16	0.17	1.25	0.17	0.12	0	0.11	0.017
18	R	0.39	0.04	0.53	0.03	1.25	0.1	1.61	0.14	0.12	0	0.07	0
19	R	0.050 <sup>d</sup>	0.13	0.30	0.11	0.93	0.15	1.79	2.1 <sup>d</sup>	0.11	0.11	0.10	0
20	N	0.32	0	0.54	0.03	1.40	0.14	1.34	0.18	0.14	0	0.12	0
21	N	lost <sup>f</sup>	lost <sup>f</sup>	0.5	0	0.8	tr <sup>e</sup>	1.0	tr <sup>e</sup>	0.02	0	0.2	0
22	R	0.46	0	0.52	0	1.28	0.13	1.16	0.14	0.14	0	0.14	0
23	R	0.415	0.03	0.535	0.04	2.035	0.189	1.365	0.16	0.085	0.01	0.175	0.02
25	R	0.25	0.043	0.48	0.037	1.38	0.12	1.3	0.11	0.10	0	0.11	0
26	R	0.40	0.06	0.44	0.04	1.35	0.13	1.24	0.12	0.12	0	0.12	0
27	R	0.41	0.03	0.44	0.04	1.18	0.11	1.30	0.13	0.11	0	0.11	0
33	R	0.50	0.08	1.00 <sup>d</sup>	0.13	1.75	0.25	1.60	0.25	0.21	0	0.24	0

\*Duplicate pairs: 1/8; 2/9; 3/10; 5/12; 6/13; 7/14.

<sup>b</sup>Determined by LC method of Foos and Warren (11). N = normal phase; R = reverse phase.<sup>c</sup>Values omitted from calculations after applying Youden's ranking test (16).<sup>d</sup>Values omitted from calculations as outlined by Dixon's test (15).<sup>e</sup>Collaborator reported trace. Trace was taken as 0.005 ng M<sub>2</sub>/mL for statistical purposes.<sup>f</sup>Sample lost.

M<sub>1</sub>; samples 4/11, no aflatoxin; samples 5/12, 17.5 µL stock M<sub>1</sub> with 18.8 µL 10-fold dilution stock M<sub>2</sub>; samples 6/13, 49.4 µL stock M<sub>1</sub> with 75.2 µL 10-fold dilution stock M<sub>2</sub>; and samples 7/14, 43.6 µL 10-fold dilution stock M<sub>1</sub>. All sample solutions and LC standard solutions were dispensed into 2 mL glass ampules (1.5 mL each), and the glass ampules were sealed. The final aflatoxin concentrations in milk were 0.078–1.31 ng M<sub>1</sub>/mL and 0.030–0.13 ng M<sub>2</sub>/mL.

### Description of Studies

Each of 14 collaborators received the following items: 1 ampule of aflatoxin M<sub>1</sub>–M<sub>2</sub> standard solution (0.50 ng M<sub>1</sub> and 0.10 ng M<sub>2</sub>/mL, in acetonitrile–benzene, 1 + 9); 1 ampule of TFA; 1 ampule of dichlorodimethylsilane (DDS); 1 ampule of practice milk-contaminating solution (0.77 ng M<sub>1</sub>, and 0.06 ng M<sub>2</sub>/mL milk); 14 coded ampules containing milk-contaminating solutions (in acetonitrile); 5 polypropylene Econo-Columns with 35 µm support disk (Bio-Rad); 17 Sep-Pak C<sub>18</sub> cartridges (Waters Associates); 20 g E. Merck silica gel 60, particle size 0.040–0.063 mm (No. 9385) containing 1% water; 5 disposable syringe tips for rubber stopper vacuum filtration; two 30 mL polypropylene syringes; and a copy of study instructions, method description, and report sheet. Collaborators were required to furnish domestic aflatoxin-free milk to prepare the artificially contaminated milk samples. For each sample, collaborators were to quantitatively pipet 1.0 mL acetonitrile solution from a sample ampule into a 100 mL graduate containing about 50 mL milk. They were to add additional milk to the 100 mL mark, pour the spiked milk into a beaker to mix, and use 20 mL of the solution for analysis as described in the method.

This procedure diluted the acetonitrile solution to prevent early elution of M<sub>1</sub> and/or M<sub>2</sub> from the C<sub>18</sub> Sep-Pak cartridges during the extraction step.

Sample solutions were prepared which would test the method at levels that have been reported in commercial milks. Each collaborator's samples were assigned a different set of computer-selected random numbers from 1 to 14. Each sample had a duplicate in the set (7 sets of blind duplicates).

### Aflatoxins M<sub>1</sub> and M<sub>2</sub> in Fluid Milk

#### Liquid Chromatographic Method

##### First Action

#### 26.B02

#### Principle

Aflatoxins M<sub>1</sub> and M<sub>2</sub> are extd from milk on C18 cartridge, eluted with ether onto silica column, eluted with CH<sub>2</sub>Cl<sub>2</sub>–alcohol, and derivatized with trifluoroacetic acid. Liq. chromatgc peaks are detected fluorometrically and compared with std-TFA derivatives.

#### 26.B03

#### Reagents

(a) *Solvents*.—Distd in glass CH<sub>3</sub>CN, CH<sub>2</sub>Cl<sub>2</sub>, and isopropyl alcohol; reagent grade alcohol, ether (0.01% EtOH preservative), hexane, MeOH, trifluoroacetic acid, and H<sub>2</sub>O (deionized, filtered thru 0.45 µm filter).

(b) *Water–acetonitrile wash soln.*—95 + 5.

(c) *Methylene chloride–alcohol elution soln.*—95 + 5.

(d) *Mobile phase*.—Prep. H<sub>2</sub>O–isopropyl alcohol–CH<sub>3</sub>CN (80 + 12 + 8). Degas in ultrasonic bath, or equiv. Alternative solv. proportions may be used to give optimum resolution (i.e., 84 + 11 + 5).

Table 2. Statistical evaluation of results for LC determination of aflatoxins M<sub>1</sub> and M<sub>2</sub> in artificially contaminated fluid milk<sup>a</sup>

Statistic	Sample 1/8	Sample 2/9	Sample 3/10	Sample 5/12	Sample 6/13	Sample 7/14
Aflatoxin M <sub>1</sub>						
Mean, ng/mL	0.656 (0.603) <sup>b</sup>	0.155 (0.143)	0.072 (0.067)	0.412 (0.415)	1.235 (1.300)	0.123 (0.116)
Std dev.	0.131	0.043	0.016	0.079	0.243	0.024
Coeff. of var., %	19.9	27.7	22.2	19.2	19.7	19.5
Theoretical, ng/mL	0.77	0.155	0.078	0.47	1.31	0.12
Recovery, %	85.2	100.0	92.3	87.7	94.3	102.5
N <sup>c</sup>	25 (20)	26 (20)	22 (18)	23 (18)	26 (20)	26 (20)
Aflatoxin M <sub>2</sub>						
Mean, ng/mL	0.062 (0.065)			0.038 (0.050)	0.128 (0.134)	
Std dev.	0.033			0.016	0.016	
Coeff. of var., %	53.2			42.1	12.5	
Theoretical, ng/mL	0.06			0.03	0.13	
Recovery, %	103.3			126.7	99.5	
N <sup>c</sup>	26 (20)			25 (20)	25 (19)	

<sup>a</sup>Calculated from values in Table 1.<sup>b</sup>Values calculated after omitting data obtained using normal phase LC (Collaborators 8, 20, 21).<sup>c</sup>N = number of values.

(e) *Aflatoxin std solns.*—Aflatoxin M<sub>1</sub> (Eureka Laboratories, Sacramento, CA 95816) and aflatoxin M<sub>2</sub> (Sigma Chemical Co.). Prep. stock solns (ca 200 µg M<sub>1</sub>/mL and 100 µg M<sub>2</sub>/mL) in CH<sub>3</sub>CN and det. concns according to 26.004–26.011, using extinction coefficients of 19 850 and 21 400 for M<sub>1</sub> and M<sub>2</sub>, resp., in CH<sub>3</sub>CN. Make working std soln contg 0.50 µg M<sub>1</sub> and 0.10 µg M<sub>2</sub>/mL in CH<sub>3</sub>CN–benzene (1 + 9) for use in prep M<sub>1</sub>-TFA derivative.

(f) *Dichlorodimethylsilane (DDS).*—5% in toluene. Add 5 mL DDS (99%) (Aldrich Chemical Co., or equiv.) to toluene and dil. to 100 mL. Store in g-s flask in cold. (Caution: DDS is a lachrymator and is flammable.)

## 26.B04

## Apparatus

(a) *Silica gel cleanup columns.*—0.8 × 4.0 cm polypropylene Econo-Column with Luer tip, 35 µm, porous polypropylene bed support disk, and 10 mL reservoir (Bio-Rad Laboratories, Cat. No. 731-1550).

(b) *Silica gel cleanup column packing and preparation.*—Dry silica gel 60, particle size 0.040–0.063 mm (E. Merck, No. 9385) in 105° oven for 1 h. Cool and add 1% H<sub>2</sub>O by wt. Shake in sealed container and equilibrate overnight before use. Assemble polypropylene column and 250 mL vac. flask fitted with 1-hole stopper as shown in Fig. 26:B1. Fill column to ca 2 mL mark with silica gel (ca 1 g). Pull gentle vac. to pack bed and add ca 1 g anhyd. Na<sub>2</sub>SO<sub>4</sub> to top of silica gel bed.

(c) *Extraction cartridges.*—C18 Sep-Pak sample prepn cartridges (Waters Associates, Inc.).

(d) *Disposable pipet tips.*—50 and 200 µL Eppendorf or equiv.

(e) *Liquid chromatograph.*—Any pulse-free or pulse-dampened liq. chromatgc system which includes pump(s), injector, and compatible recorder.

(f) *Fluorescence detector.*—Any fluorescence detector capable of providing 365 nm excitation and >400 nm emission wavelengths and sensitivity of 50–100% full-scale response for 1 ng M<sub>1</sub>-TFA derivative (e.g., Kratos-Schoeffel FS 970).

(g) *LC analytical column.*—Any 0.4 × 25 cm column contg spherical 5 µm particle size C18 bonded silica gel (e.g., DuPont ODS, Spherisorb 5 ODS II).

(h) *Vacuum regulator.*—Any com. or custom device capable of creating and controlling partial and full vac. with side arm vac. flask.

(i) *Silylated vials for aflatoxin std solns.*—Fill 1 or 1½ dram glass vials nearly full with 5% DDS and heat ca 40 min at 45–55°. Discard soln, and rinse vials 3 times with toluene and then 3 times with MeOH. Heat vials in oven at 75° for 20–30 min to evap. MeOH. Cap vials (with Teflon liners) and store for aflatoxin std solns.

## 26.B05

## Extraction

Attach inlet (longer) stem of C18 cartridge to Luer tip of 30–50 mL syringe. Assemble syringe, cartridge, and vac. flask as shown in Fig. 26:B1. Adjust vac. to pull solvs thru cartridge in fast *dropwise* manner (ca 5 mm Hg). Prime cartridge by adding 5 mL MeOH, then 5 mL H<sub>2</sub>O (do not pull cartridge dry; leave small excess H<sub>2</sub>O in stem). Discontinue vac. and remove cartridge-syringe assembly from stopper to prevent loss of prime.

Warm sample to room temp. Gently invert sample ≥10 times to evenly distribute cream in nonhomogenized samples. Transfer 20 mL milk to graduate contg 20 mL hot (ca 80°) H<sub>2</sub>O. (If necessary, more hot H<sub>2</sub>O may be used to thin milk soln.)

Replace cartridge-syringe assembly in stopper. Pour entire 40 mL warm, dild sample into syringe and gently pull sample thru cartridge at flow rate ca 30 mL/min (very fast drops). *Caution:* Too fast a flow will not allow sufficient time for aflatoxin to adsorb, resulting in low recoveries. Add 10 mL H<sub>2</sub>O–CH<sub>3</sub>CN wash soln to syringe and pull thru. Plug syringe barrel with stopper and pull hard vac. on cartridge for ca 30 s to remove as much wash soln as possible from packing. Remove cartridge and dry inside of both stems with cotton swab or tissue paper to eliminate any remaining wash soln. Reprime cartridge by adding 150 µL CH<sub>3</sub>CN to inlet bed support disk and let solv. soak into packing for 30 s. Attach cartridge to dry glass or plastic 10 mL Luer tip syringe, retaining same stem as inlet.

Insert silica gel cleanup column into 250 mL vac. flask fitted with 1-hole rubber stopper (Fig. 26:B1). Wash column with 5 mL ether. Add 7 mL ether to syringe-cartridge positioned above silica gel cleanup column. With plunger, slowly force ether thru cartridge (in portions), collecting eluate in column reservoir. Pull ether slowly thru silica cleanup column, using vac. to maintain flow rate ca 10 mL/min (fast drops). Rinse silica column with 2 mL addnl ether, continuing to use vac. Discard ether.

Remove column and stopper from flask and place 16 × 125 mm collection tube in flask to catch eluate from column. Add 7 mL




FIG. 26:B1—Diagram of apparatus for extraction and cleanup of milk extracts

**Table 3. Collaborative results for determination of aflatoxins M<sub>1</sub> and M<sub>2</sub> (ng/mL) in uncontaminated fluid milk (blind duplicates)<sup>a</sup>**

Coll.	Method <sup>b</sup>	Sample 4		Sample 11	
		M <sub>1</sub>	M <sub>2</sub>	M <sub>1</sub>	M <sub>2</sub>
1 <sup>c</sup>	R	0.10	<0.01	0.02	0
4	R	0	0	0	0
6	R	0	0	0	0
8	N	0	0	0	0
18	R	0	0	0.03	0
19	R	0	0	0	0
20	N	0	0	0.02	0
21	N	0.08	0	tr <sup>d</sup>	0
22	R	0	0	0	0
23	R	0	0	0	0
25	R	0	0	0	0
26	R	tr	0	0	0
27	R	0	0	0	0
33	R	0	0	0	0

Method	Sample 4		Sample 11		Other samples <sup>e</sup>
	M <sub>1</sub>	M <sub>2</sub>	M <sub>1</sub>	M <sub>2</sub>	M <sub>2</sub>
Total observations:	13	13	13	13	78
Positive observations:					
Reverse phase	1	0	1	0	11
Normal phase	1	0	2	0	3
Negative observations:					
Reverse phase	9	10	9	10	49
Normal phase	2	3	1	3	15

<sup>a</sup>Determined by LC method of Foos and Warren (11).<sup>b</sup>N = normal phase; R = reverse phase.<sup>c</sup>Values omitted from calculations after applying Youden's ranking test (16).<sup>d</sup>Collaborator reported trace.<sup>e</sup>Samples 2/9, 3/10, and 7/14.

elution solv. (CH<sub>2</sub>Cl<sub>2</sub>-alcohol) to column reservoir. Pull solv. thru column with vac. at ca 10 mL/min flow rate, collecting eluate in tube.

Discontinue vac. and remove collection tube from assembly. Evap. eluate just to dryness under N stream, using heat to keep collection tube near room temp. or under vac. at <35°.

Transfer residue to 1 dram vial with CH<sub>2</sub>Cl<sub>2</sub> and evap. to dryness under N on steam bath or in heating block <50°. (Do not overheat dry sample.) Save for derivative prepn.

**26.B06****Liquid Chromatography**

Prep. derivative of sample exts by adding 200 µL hexane and 200 µL trifluoroacetic acid to dry residue in vial. Shake on vortex mixer ca 5–10 s. Let mixt. sit 10 min at 40°, in heating block or bath; then evap. to dryness under N on steam bath or heating block (<50°). Add 2 mL H<sub>2</sub>O-CH<sub>3</sub>CN (75 + 25) to vial to dissolve residue and shake well in vortex mixer for LC analysis. For derivative of std M<sub>1</sub>, add 200 µL hexane and 50 µL trifluoroacetic acid to silylated vial and mix. Add 50 µL M<sub>1</sub>-M<sub>2</sub> working std soln *directly* into hexane-TFA mixt. and shake on vortex mixer 5–10 s. Treat as described for sample derivative. Stabilize instrument and detector for suitable period at flow rate of 1.0 mL/min with H<sub>2</sub>O-isopropanol-CH<sub>3</sub>CN (80 + 12 + 8). Adjust detector attenuator so that 50–100 µL injection of std (0.625–1.25 ng M<sub>1</sub>, 0.125–0.25 ng M<sub>2</sub>) gives 50–75% full-scale recorder pen deflection for aflatoxin M<sub>1</sub>. Inject LC std 2–3 times until peak hts are const. Prep. std curve from either peak hts or peak areas to ensure linear relationship. Inject sample exts (typically 50–100 µL) with std injections interspersed to ensure accurate quantitation. Retention times of M<sub>1</sub> (as TFA derivative) and M<sub>2</sub> are ca 4–5 min and ca 7 min, resp.

Calc. aflatoxin concn:

$$\text{ppb (M}_1 \text{ or M}_2\text{)} = (H \times C' \times VI' \times V)/(H' \times VI \times W)$$

where *H* and *H'* = peak ht of sample and std, resp.; *C'* = concn of std (ng/µL); *VI'* and *VI* = vol. injected of std and sample, resp.;

*V* = final total sample vol. (µL); and *W* = vol. of milk represented by final ext (typically 20 mL). Sep. calc. concn for M<sub>1</sub> and M<sub>2</sub>.

**Normal Phase LC Procedure**

The normal phase LC procedure, which was not adopted, was performed as follows:

**Reagents**

(a) *Solvents*.—Distilled-in-glass acetonitrile, chloroform, and methylene chloride; reagent grade diethyl ether (0.01% ethanol preservative), ethanol, and methanol.

(b) *Water-acetonitrile wash solution*.—95 + 5.

(c) *Methylene chloride-ethanol elution solution*.—95 + 5.

(d) *Mobile phase*.—Add 22.5 mL ethanol (reduce by amount in CHCl<sub>3</sub> as preservative) to 1 L CHCl<sub>3</sub> solution (225 mL water-saturated CHCl<sub>3</sub> + 775 mL CHCl<sub>3</sub> [ethanol-preserved]). Stir and degas in ultrasonic bath or equivalent.

(e) *Aflatoxin standard solutions*.—Aflatoxin M<sub>1</sub> (Eureka Laboratories, Sacramento, CA 95816) and aflatoxin M<sub>2</sub> (Sigma Chemical Co., St. Louis, MO 63176). Prepare separate stock solutions of each in acetonitrile-benzene (1 + 9) to concentration of 1.0 µg/mL. Store in freezer. Prepare standard daily by diluting aliquots of stock solutions with fresh mobile phase to obtain 1 solution that contains 0.01 µg M<sub>1</sub>/mL and 0.004 µg M<sub>2</sub>/mL.

(f) *Dichlorodimethylsilane (DDS)*.—See 26.B03(f).

**Apparatus**

(a) *Silica gel cleanup columns*.—See 26.B04(a).

(b) *Silica gel cleanup column packing and preparation*.—See 26.B04(b).

(c) *Extraction cartridges*.—See 26.B04(c).

(d) *Disposable pipet tips*.—See 26.B04(d).

(e) *Liquid chromatograph*.—See 26.B04(e).

(f) *Fluorescence detector*.—Any fluorescence detector with silica gel-packed cell and filters to provide 365 nm excitation and >400 nm emission wavelengths and sensitivity of 50–100% full-scale response for 1 ng M<sub>1</sub> (e.g., Varian Fluorichrom).

(g) *LC analytical column*.—Any 0.4 × 25 cm column containing spherical-shape 5 µm particle size silica gel (e.g., DuPont Zorbax Sil, Varian Micro-Pak SI-5).

(h) *Vacuum regulator*.—See 26.B04(h).

(i) *Silylated vials for aflatoxin standard solutions*.—See 26.B04(i).

**Extraction**

See 26.B05 through next-to-last paragraph. Then dissolve eluate residue in 1000 µL fresh LC mobile phase and mix well.

**Liquid Chromatography**

Stabilize instrument and detector for suitable period at flow rate of 1.0 mL/min with CHCl<sub>3</sub>-ethanol mobile phase. Adjust detector attenuator so that 100 µL injection of standard (1.0 ng M<sub>1</sub>, 0.4 ng M<sub>2</sub>) gives 50–75% full-scale recorder pen deflection for aflatoxin M<sub>1</sub>. Inject 100 µL LC standard 2–3 times or until peak heights are constant. Prepare standard curve from either peak heights or peak areas to ensure linear relationship. Inject sample extracts (typically 100 µL) with standard injections interspersed to ensure accurate quantitation. Retention times of M<sub>1</sub> and M<sub>2</sub> are 8–10 min and 9–12 min, respectively, depending on length of column.

Calculate aflatoxin concentration as shown in 26.B06.

Table 4. Analysis of variance for combined samples<sup>a</sup>

Source of variation	Aflatoxin M <sub>1</sub>		Aflatoxin M <sub>2</sub>		Expected mean square
	df <sup>b</sup>	Mean square	df	Mean square	
Laboratories (L)	9	0.19073***	9	0.26283***	S <sub>L</sub> <sup>2</sup> + 2S <sub>Ls</sub> <sup>2</sup> + 12S <sub>L</sub> <sup>2</sup> (M <sub>1</sub> ) or 6S <sub>L</sub> <sup>2</sup> (M <sub>2</sub> )
Samples (S)	5	4.78924***	2	0.98834***	
L × S	45	0.02105 ns	18	0.03335***	S <sub>S</sub> <sup>2</sup> + 2S <sub>Ls</sub> <sup>2</sup>
Duplicates	60	0.01145	30	0.00867	S <sub>S</sub> <sup>2</sup>

Component	Precision Parameters, % <sup>c</sup>	
	M <sub>1</sub>	M <sub>2</sub>
Repeatability (S <sub>o</sub> )	27.9	23.9
Lab-sample interaction (S <sub>Ls</sub> )	17.3	29.1
Among labs (S <sub>L</sub> )	31.5	56.9
Reproducibility <sup>d</sup> (S <sub>x</sub> )	44.5	64.7

\*\*\* = Significant at 0.001 level; ns = not significant

<sup>a</sup>Analysis of variance calculated by log transformation according to Snedecor and Cochran (18). Values are based on logarithms. Normal phase data (Cols. 8, 20, and 21) not included.

<sup>b</sup>df = degrees of freedom

<sup>c</sup>Calculated from antilogarithm formula:  $(10 \sqrt{S^2} - 1) \times 100$ ; for example,  $S_{Ls} = (10 \sqrt{\frac{0.02105 - 0.01145}{2}} - 1) \times 100 = 17.3\%$

<sup>d</sup> $S_x^2 = S_o^2 + S_L^2$

### Results and Discussion

Individual values were omitted from calculations according to Dixon's test for outliers at the 0.05 level (15) using either upper or lower 1-tail test (never both on same sample). For statistical calculations based on log transformations, the retained value of the sample pair was substituted for the outlier to maintain balance in the analysis of variance. The "lost" sample (Collaborator 21, sample 5/12) was treated similarly. For aflatoxin M<sub>2</sub> values reported as "trace," a value of 0.005 ng/mL was substituted. The values for Collaborator 1 were not included in the calculations because the composite data exceeded the lower limit of Youden's ranking test (16). The decision to omit these data was made because 11 of the 12 sample values were either the lowest value (7 samples) or the second lowest (4 samples) reported by the collaborators. This was an abnormally high percentage of low values.

The results reported for aflatoxins M<sub>1</sub> and M<sub>2</sub> in artificially contaminated fluid milk are shown in Table 1, and the statistical summaries are given in Table 2. Only 3 collaborators used normal phase LC for their final determination; consequently, separate statistical comparisons between normal and reverse phase LC data were not made. An analysis of variance was done with the reverse phase data only because normal phase data are insufficient to give meaningful results.

The statistical means for the 6 duplicate sets of artificially contaminated samples are shown in Table 2. Aflatoxin M<sub>1</sub> concentrations ranged from 0.072 to 1.235 ng/mL (ppb), while aflatoxin M<sub>2</sub> concentrations varied from 0.038 to 0.128 ng/mL. Levels for both toxins are similar to those reported in contaminated commercial milk. No false negative values were reported for the lowest M<sub>1</sub> concentration (sample 3/10, 0.072 ng/mL); therefore, a minimum detection limit of 70 ppt is realized without modifications of sample extract volumes. Lower detection limits (10–50 ppt) have been obtained (at NRRC) by reducing the final extract volume from 2000 µL to 500–1000 µL.

Aflatoxin M<sub>1</sub> recoveries were excellent (Table 2) and ranged from 85.2 to 102.5%. The average, 93.7%, compares favorably with that of the previous collaborative study (1) (method I, 91% recovery). The standard deviations for M<sub>1</sub> samples are considerably less in this study as are the coefficients of variation (Table 2) (19.2–27.7%, av. 21.4%) compared with those of method I (47%) (1).

Five false positive values, including 2 "trace" values, were reported in this current study for 26 uncontaminated M<sub>1</sub> samples (Table 3). This is comparable to the 2 previous AOAC aflatoxin M<sub>1</sub> collaborative studies (1, 17); however, a lesser number of false positives was expected. Two collaborators (Collaborators 18 and 23) necessarily used contaminated commercial milk and subtracted an "average" background value from their sample concentrations. One reported a false positive (Collaborator 18, sample 4/11, Table 3). Low-level contamination in commercial milk (<0.09 ng M<sub>1</sub>/mL) had been determined by several scientists in 1984 (private communications). This problem may have contributed to the larger incidence of normal phase false positives (3 of 6).

Recoveries of aflatoxin M<sub>2</sub> from the artificially contaminated milk were elevated slightly: 99.5, 103.3, and 126.7% (Table 2). The highest recovery (samples 5/12, 126%) was for the lowest level samples, 0.03 ng/mL. Peak heights or areas for aflatoxin M<sub>2</sub> are not large for low concentrations; therefore, a mean of 0.038 ng/mL for a theoretical value of 0.03 ng/mL is satisfactory and acceptable. Aflatoxin M<sub>2</sub> coefficients of variations were twice those for M<sub>1</sub> except for samples 6/13 which were lower (12.5%). This reflects the increased level of M<sub>2</sub> in samples 6/13 (0.128 ng/mL). Since no other study has included aflatoxin M<sub>2</sub>, comparisons are not possible. Very little data have been published for aflatoxin M<sub>2</sub> in commercial samples; so, either M<sub>2</sub> is not a contamination problem or scientists are not looking for or reporting it. Of the 76 positive aflatoxin M<sub>2</sub> samples, there were only 6 false negatives (<10%) (Table 1). Four were from the lowest level samples (samples 5/12, 0.03 ng/mL). There were 13 false positives for 78 recorded values (Tables 1 and 3). Eight were found on only one of the possible two in each duplicate set of negative samples. This result indicates that care needs to be exercised when identifying aflatoxin M<sub>2</sub> peaks in milk samples by LC. There were no false positive aflatoxin M<sub>2</sub> values for the uncontaminated milk samples 4/11 (Table 3). Evidently, collaborators realized that study samples without aflatoxin M<sub>1</sub> were not likely to have aflatoxin M<sub>2</sub> either.

The analyses of variance and variance components after combining data for all artificially contaminated samples from reverse phase data are presented in Table 4. A better understanding of the total variance and the component sources of variation is obtained by use of the entire composite of samples. A log transformation was used to compute the analysis

of variance as recommended by Snedecor and Cochran (18) for situations where standard deviations (Table 2) are proportional to means rather than constant. This relation was observed in this study for although the CV values were relatively constant, sample means varied by a factor of about 20. The symbols representing the precision parameters are the ones recommended by the Committee on Collaborative Interlaboratory Studies (19).

The within-laboratory precision parameter,  $S_o$ , is the repeatability while the among-laboratory precision parameter,  $S_x$ , is the reproducibility. The latter includes both the within-laboratory variance,  $S_o$ , and the among-laboratory variance,  $S_L$ . For aflatoxin  $M_1$ , the repeatability ( $S_o$ ) was 27.9% which compares favorably with a previous study, 25% (1). The reproducibility, 44.5%, was less than calculated for Method I (47%) (1). It is noteworthy that the lab-sample interaction values are very low for both  $M_1$  and  $M_2$ , even though collaborators provided their own fluid milk. The repeatability for aflatoxin  $M_2$  was 23.9% whereas the reproducibility was 64.7%. The latter value is higher than for  $M_1$ ; however, there is no other  $M_2$  study for comparison.

Generally, the collaborators commented positively about the method. They liked its rapidity and the lack of chromatographic interferences. Only one collaborator (No. 19) did not feel the method was rapid. Four collaborators (Nos. 18, 19, 22, and 25) detected incomplete derivative reactions in the reverse phase LC procedure. Subsequently, these collaborators tried heating the reaction mixture and found no unreacted aflatoxin  $M_1$ . Collaborator 22 suggested that the method should incorporate a statement to store the dry film extracts in the freezer until the LC step. Collaborator 25 commented that silylating the vials helped prevent degradation of aflatoxin standards and was in favor of this technique.

### Recommendation

The conclusion from the evaluation of this collaborative study indicates that the Foos-Warren method for the rapid determination of aflatoxins  $M_1$  and  $M_2$  in fluid milk by reverse phase LC (11) should be adopted as official first action. The Associate Referee on aflatoxin M recommends that action. Although the data for normal phase LC are satisfactory, they are insufficient to give a meaningful analysis of variance; therefore, normal phase LC cannot be recommended for adoption.

### Acknowledgments

The authors thank W. Bailey for statistical processing of the data; M. Sloan, formerly of Waters Associates, Milford, MA, for providing the  $C_{18}$  Sep-Pak cartridges; and the following collaborators for their cooperation:

M. Billotte, United Dairymen Association of Arizona, Tempe, AZ

J. Ferguson-Foos and J. D. Warren, Florida Dept of Agriculture and Consumer Services, Tallahassee, FL

J. I. Greer, U.S. Dept of Agriculture, Northern Regional Research Center, Peoria, IL

K. D. Johnson, General Mills, Inc., Minneapolis, MN  
 J. I. Pitt and N. F. Tobin, Commonwealth Scientific and Industrial Research Organization, North Ryde, Australia  
 J. L. Richard and R. Lyon, U.S. Dept of Agriculture, National Animal Disease Center, Ames, IA  
 B. Roberts, Central Veterinary Laboratory, Ministry of Agriculture, Fisheries, and Food, New Haw, Weybridge, Surrey, UK  
 Y. Saito, National Institute of Hygienic Sciences, Division of Foods, Tokyo, Japan  
 G. Sutton, South Carolina Dept of Health and Environmental Control, Columbia, SC  
 M. Trucksess, Food and Drug Administration, Washington, DC  
 H. P. van Egmond and E. A. Sizoo, National Institute of Public Health, Bilthoven, The Netherlands  
 N. Webb and M. E. Rhodes, Florida Dept of Agriculture and Consumer Services, Tallahassee, FL  
 D. Wilson and R. Beaver, University of Georgia, Dept of Plant Pathology, Coastal Plains Experiment Station, Tifton, GA  
 A. E. Yousef, G. Rusul, and E. H. Marth, Dept of Food Science, University of Wisconsin, Madison, WI

### REFERENCES

- (1) Stubblefield, R. D., van Egmond, H. P., Paulsch, W. E., & Schuller, P. L. (1980) *J. Assoc. Off. Anal. Chem.* **63**, 907-921
- (2) Stubblefield, R. D. (1979) *J. Am. Oil Chem. Soc.* **56**, 800-802
- (3) van Egmond, H. P., Paulsch, W. E., & Schuller, P. L. (1978) *J. Assoc. Off. Anal. Chem.* **61**, 809-812
- (4) "Changes in Methods" (1980) *J. Assoc. Off. Anal. Chem.* **63**, 309
- (5) Winterlin, W., Hall, G., & Hsieh, D. P. H. (1979) *Anal. Chem.* **51**, 1873-1874
- (6) Beebe, R. M., & Takahashi, D. M. (1980) *J. Agric. Food Chem.* **28**, 481-482
- (7) Fremy, J. M., & Boursier, B. (1981) *J. Chromatogr.* **219**, 156-161
- (8) Tuinstra, L. G. M. T., & Haasnoot, W. (1982) *Fresenius Z. Anal. Chem.* **312**, 622-623
- (9) Chang, H. L., & DeVries, J. W. (1983) *J. Assoc. Off. Anal. Chem.* **66**, 913-917
- (10) Cohen, H., LaPointe, M., & Fremy, J. M. (1984) *J. Assoc. Off. Anal. Chem.* **67**, 49-51
- (11) Ferguson-Foos, J., & Warren, J. D. (1984) *J. Assoc. Off. Anal. Chem.* **67**, 1111-1114
- (12) *Official Methods of Analysis* (1980) 13th Ed., AOAC, Arlington, VA
- (13) Stubblefield, R. D., Shotwell, O. L., & Shannon, G. M. (1968) *J. Am. Oil Chem. Soc.* **45**, 686-688
- (14) Stubblefield, R. D., Shannon, G. M., & Shotwell, O. L. (1971) *J. Am. Oil Chem. Soc.* **47**, 389-390
- (15) Dixon, W. J., & Massey, F. J. (1957) *Introduction to Statistical Analysis*, 2nd Ed., McGraw-Hill Book Co., New York, NY, p. 277
- (16) Youden, W. J., & Steiner, E. H. (1975) *Statistical Manual of the AOAC*, AOAC, Arlington, VA
- (17) Stubblefield, R. D., & Shannon, G. M. (1974) *J. Assoc. Off. Anal. Chem.* **57**, 852-857
- (18) Snedecor, G. W., & Cochran, W. C. (1980) *Statistical Methods*, 7th Ed., Iowa State University Press, Ames, IA, pp. 291-292
- (19) Horwitz, W. (1984) *J. Assoc. Off. Anal. Chem.* **67**, 432