# INADEQUACY OF SELECTIVE PLATING MEDIA IN FIELD DETECTION OF SALMONELLA

N. A. COX<sup>1</sup> and M. E. BERRANG

USDA, ARS, Russell Research Center, Athens, GA 30604-5677

Phone: (706) 546-3484

FAX: (706) 546-3772

E-mail: ncox@ars.usda.gov

Primary Audience: Quality Control Personnel, Researchers, Laboratory Personnel

# SUMMARY

Salmonella is an important human enteropathogen associated with poultry. Poultry is routinely analyzed for Salmonella by industry, regulatory, and research laboratories. Effectiveness of selective plating media is critical to Salmonella laboratory analysis. Many media routinely used in the food microbiology laboratory were developed for all Gramnegative human enteric pathogens encountered in a clinical laboratory, rather than for Salmonella alone. The presence of extraneous, non-Salmonella bacteria make the recognition of Salmonella on these plates very difficult. This, more than any other factor, contributes to poor performance of a medium. Primarily because of the poor performance of many selective agars in detecting Salmonella from non-clinical samples, food microbiology laboratories usually use two or more plating media to reduce false negative results.

Key words: Bacteria, carcass rinse, media

2000 J. Appl. Poultry Res. 9:403-406

# **DESCRIPTION OF PROBLEM**

The laboratory isolation of Salmonella began approximately a century ago. Since that time there have been numerous reports in the scientific literature comparing the efficacy of various media and incubation conditions. An important part of the Salmonella laboratory procedure is the selective plating agar. Many of the currently used media were developed for the clinical microbiologist, and clinical samples have been used for the initial evaluations of these media [1]. Although these plating media were not initially developed for isolating Salmonella from environmental

sources (and more specifically, poultry sources), several published reports have evaluated selective plating media efficacy with poultry samples [2, 3, 4, 5]. However, the conflicting results of these past studies render them of limited value to the Salmonella researcher. It is rare to find any two papers that recommend the same plating medium on the basis of efficacy. One report may indicate that a certain plating medium is best, and another will report that the same plating medium was the poorest performer of those tested. Therefore, the objective of this paper is not to compare media, but rather to demonstrate the inability of many commonly

<sup>&</sup>lt;sup>1</sup> To whom correspondence should be addressed

used plating media to adequately detect naturally occurring and artificially inoculated Salmonella from poultry carcass rinses. Food and poultry microbiologists need to be aware of the inherent error and limitations of their laboratory analyses of various food products and the primary reasons for these limitations.

## MATERIALS AND METHODS

# SAMPLES WITH NATURALLY OCCURRING SALMONELLA

A sample of 24 broiler carcasses was collected directly from the end of the chill tank of a commercial processing plant. Each carcass was placed in a separate sterile plastic bag, 100 mL of sterile distilled water was added [6], and the contents were vigorously shaken for 1 min in a 1-ft arc [4]. The resulting rinse fluid was then poured into sterile specimen cups, packed in crushed ice, and transported to the laboratory. Concentrated (10X) Universal Pre-enrichment broth (UP; Difco, Sparks, MD) was added to the rinse fluid, vielding a single-strength UP. The UP broth tubes were incubated at 35°C for 24 hr, then 0.1 mL of UP was aseptically transferred to TT Haina broth (Difco, Sparks, MD), which was incubated at 42°C for 24 hr. Next, one of each of the selective plates used in this study (Table 1) was streaked to produce isolated colony-forming units. All plates were made fresh and stored at 4°C until use. After 24 hr incubation of the plates at 37°C, three different typical Salmonella colony-forming units were picked from each type of media and used to inoculate Triple Sugar Iron (TSI) and Lysine Iron (LIA) agar (Difco, Sparks, MD) slants. Cultures giving typical reactions in the TSI and LIA tubes were confirmed to be Salmonella using Salmonella Poly O and Poly H serology (Difco, Sparks, MD).

#### SAMPLES WITH ARTIFICIALLY INOCU-LATED SALMONELLA

Approximately 10<sup>6</sup> cells of Salmonella typhimurium were inoculated into 100 mL of chicken rinse fluid (which provided a mixed background microflora), thoroughly agitated, and then streaked onto four different selective agar plates (BGS, MLIA, RP, and XLD) to produce isolated colony-forming units. After streaking, the plates were incubated at 37°C for 24 hr. Three separate colony-forming units

TABLE 1. Ability of seven different selective plating media to detect the presence of naturally occurring salmonellae from processed broiler carcass rinse samples

| SELECTIVE<br>PLATING MEDIA<br>USED             | SALMONELLAE-<br>POSITIVE<br>ISOLATIONS <sup>A</sup> |    |
|--|---|----|
|  | #   | %  |
| BG sulfa (BGS) <sup>B</sup>                    | 17  | 89 |
| Bismuth sulfite (BS) <sup>B</sup>              | 1   | 5  |
| Hektoen enteric (HE) <sup>B</sup>              | 11  | 58 |
| Modified lysine iron (MLIA)                    | 17  | 89 |
| Rambach plus (RP) <sup>C</sup>                 | 16  | 84 |
| Xylose lysine desoxycholate (XLD) <sup>B</sup> | 6   | 32 |
| Xylose lysine T4<br>(XLT4) <sup>C</sup>        | 11  | 58 |
| I .  |   |    |

 $<sup>^{</sup>A}N = 24.$ 

typical of Salmonella were selected from each plate, biochemically characterized, and serologically confirmed to be the inoculated Salmonella typhimurium.

# RESULTS AND DISCUSSION

The ability of the seven different selective plating media to detect naturally occurring salmonellae from processed broiler carcass rinses is presented in Table 1. Nineteen out of 24 (79%) broiler carcass rinse samples were confirmed to be Salmonella positive with one or more of the plating agars (data not shown in table). BGS, MLIA, and RP came the closest to allowing detection of these 19 naturally contaminated Salmonella-positive samples. Any two of these three plating media would have detected 18 of the 19 positive samples. Two other plating media (HE, XLT4) displayed intermediate performance, allowing detection of these naturally occurring Salmonella, while XLD and BS performed very poorly. Other studies [7, 8, 9] also reported that HE and XLD were less effective than other media tested, although a recent survey found that both of these plating media are commonly used [9]. All seven plating agars were streaked from the same selective enrichment tubes, therefore each plate had

<sup>&</sup>lt;sup>B</sup>Difco Laboratories, Division of Becton Dickinson and Co., Sparks, MD.

<sup>&</sup>lt;sup>C</sup>Unipath Ltd., Basingstoke, UK.

approximately the same number of Salmonella cells deposited on the agar surface. In spite of that, only one rinse sample was found to be Salmonella positive with all seven of the media.

Culturing of environmental samples for the presence of Salmonella poses several challenges to the microbiologist. To begin with, the number of viable Salmonella cells present in the sample taken may be very low. Worse yet, these cells may be in an injured or altered condition. But the most critical challenge of all is that these Salmonella cells in a sample such as chicken rinse will be in the midst of millions of other non-Salmonella (extraneous) bacteria. Detecting their presence thus can be very difficult. Most selective media were not developed for Salmonella alone, but rather for all Gram-negative human enteric pathogens encountered in a clinical laboratory.

A plating medium may perform adequately in a clinical laboratory, yet not detect Salmonella from environmental samples. To further dramatize how Salmonella cells or colony-forming units can be present on a plating medium and not be recognized, chicken rinse fluid samples (100 mL) were inoculated with 1 million Salmonella cells and then streaked onto four of the plating media (BGS, MLIA, RP, or XLD). The results are presented in Table 2. Even though hundreds of Salmonella cells (approximately 300) were deposited on each of these media, only RP detected the inoculated Salmonella in all eight spiked rinse samples.

Since large numbers of samples prohibit the microbiologist from selecting more than

TABLE 2. Recovery of artificially inoculated Salmonella (106 cells) from chicken rinse fluid

| CENTRALE (100 cons) Horizontal                 |                                       |  |
|--|---------------------------------------|--|
| SELECTIVE PLATING<br>MEDIA USED                | NUMBER OF<br>SALMONELLA<br>ISOLATIONS |  |
| Rambach plus (RP) <sup>A</sup>                 | 8 <sup>B</sup>                        |  |
| Modified lysine iron (MLIA) <sup>A</sup>       | 6                                     |  |
| BG sulfa (BGS) <sup>C</sup>                    | 3                                     |  |
| Xylose lysine desoxycholate (XLD) <sup>C</sup> | 2                                     |  |
|  |                                       |  |

AUnipath Ltd., Basingstoke, UK.

2 or 3 CFU/plate, Salmonella colonies must not only be present, but present in a recognizable form. Many of the media rely on color patterns based on a pH indicator, usually a fermentable carbohydrate such as lactose. The non-Salmonella or extraneous bacteria such as coliforms or Escherichia coli will ferment the lactose and produce a color determined by the indicator. Salmonella does not normally ferment lactose, so its colonies' appearance will be considerably different on the media. However, abundance of these extraneous microorganisms can color an entire section of the plate surface, thereby causing the outnumbered Salmonella to go unnoticed. Mainly because of the poor performance of many selective agars in detecting Salmonella from environmental samples, microbiology laboratories usually use two or more plating media to minimize false negative

# CONCLUSIONS AND APPLICATIONS

- 1. Salmonella may go undetected in environmental (poultry) samples because of the presence of extraneous or background types of bacteria.
- 2. The inherent fallibility in laboratory methods must be considered when analyzing data resulting from Salmonella analysis.
- 3. Even when two laboratories are using the identical procedures, variation in the results should be expected.
- 4. Although none of the plating media tested in this study performed perfectly, some of these media do perform better than others.

<sup>&</sup>lt;sup>B</sup>N = 8 spiked carcass rinse samples.

<sup>&</sup>lt;sup>C</sup>Difco Laboratories, Division of Becton Dickinson and Co., Sparks, MD.

## REFERENCES AND NOTES

- 1. Dusch, H. and M. Altwegg, 1995. Evaluation of five new plating media for isolation of Salmonella species. J. Clinical Microbiology 33:802–804.
- 2. Moats, W.A., 1978. Comparison of four agar plating media with and without added novobiocin for the isolation of salmonellae from beef and deboned poultry meat. Appl. and Environ. Microbiology 36:747-751.
- 3. Cox, N.A., B.H. Davis, J.H. Kendall, A.B. Watts, and A.R. Colmer, 1972. Salmonella in the laying hen. III. A comparison of various enrichment broths and plating media for the isolation of Salmonella from poultry feces and poultry food products. Poultry Sci. 51:1312-1316.
- 4. Cox, N.A., 1988. <u>Salmonella</u> methodology update. Poultry Sci. 67:921-927.
- 5. Waltman, W.D., A.M. Horne, and C. Dirkle, 1995. Comparative analysis of media and methods for isolating Salmonella from poultry environmental samples. Pages 1–8 in: Annual Meeting of the United States Anim. Health Assn., Reno, NV.

- 6. Cox, N.A., J.E. Thomson, and J.S. Bailey, 1981. Sampling of broiler carcasses for Salmonella with low volume water rinse. Poultry Sci. 60:768-770.
- 7. Komatsu, K.K. and L. Restraino, 1981. Determination of the effectiveness of novobiocin added to two plating media for the isolation of <u>Salmonella</u> from fresh meat products. J. Food Safety 3:183–192.
- 8. Miller, R.G., C.R. Tate, E.T. Mallinson, and J.A. Scherrer, 1991. Xylose-lysine-tergitol 4. An imposed selective agar medium for the isolation of <u>Salmonella</u>. Poultry Sci. 70:2429-2432.
- 9. Waltman, W.D. and E.T Mallinson, 1995. Isolation of <u>Salmonella</u> from poultry tissue and environmental samples: A nationwide survey. Avian Dis. 39:45-54.

#### **ACKNOWLEDGEMENT**

The authors acknowledge excellent technical assistance by Ms. Debbie Posey.